



STaR

IHE Integration Statement

Product Version: **1.86**
Revision: **3**
Date : **September 2009**

ETIAM Head Office in Europe:
ETIAM Technopole Atalante-Champeaux
2, rue du Pierre-Joseph Colin
35000 Rennes
France

Phone: +33.(0)2.99.14.33.88
Fax: +33.(0)2.99.14.33.80

E-mail: info@etiam.com
www.etiam.com

ETIAM Head Office in America:
ETIAM Corp.
185 Alewife Brook Parkway, Suite 410
Cambridge, MA 02138
USA

Toll Free: +1-(877) 384-2662 (USA & Canada)
Phone/Fax: +1-(617) 395-5809

IHE Integration Statement			
Vendor	Product Name	Version	Date
ETIAM	STaR	1.86	September 2009
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Consistent Time	Time Client	None	
Scheduled Workflow	Patient Demographics Consumer	None	
	Performed Procedure Step Manager	None	
Patient Information Reconciliation	Image Manager/Image Archive	None	
	Performed Procedure Step Manager	None	
Access to Radiology Information	Image Manager/Image Archive	None	
	Report Repository	None	
Key Image Note	Image Manager/Image Archive	None	
Simple Image and Numeric Report	Report Repository	None	
Consistent Presentation of Images	Image Manager/Image Archive	None	
Reporting Workflow	Image Manager/Image Archive	None	
Evidence Documents	Image Manager/Image Archive	None	
Cross-enterprise Document Sharing for Imaging	Imaging Document Source	None	
Retrieve Information Display	Information Source	None	
Internet address for vendor's IHE information: http://www.etiam.com/ihe			
Links to Standards Conformance Statements for the Implementation			
DICOM	http://www.etiam.com/ihe		
HL7			
Links to General Information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.org	In Japan: http://www.jira-net.or.jp/ihe-j	